

Experiment 7


Complexometric titration using **EDTA** **(Chelatometry)**

In Complexometric titration

- ** The formation of a colored complex is used to indicate the end point of a titration.
- ** It's particularly useful for the determination of a mixture of different metal ions in solution.
- ** A simple ion is transformed into a complex ion and the equivalent point is determined by using metal indicators

What is EDTA ?

- Ethylenediaminetetraacetic acid



✱

What is EDTA ?

- ** The most titrant used.
- ** It has four carboxyl groups and two amine groups that can act as electron pair donors.
- ** EDTA is tetraprotic and hexadentate ligand
- ** EDTA forms stable, water soluble 1:1 ratio complex with metal.
- ** EDTA can be obtained as H_4Y or H_2Y^{2-}
- ** H_2Y^{2-} can form only 4 coordinate covalent bonds.
- ** EDTA is extensively used to standardize metal (2+) cations



- as stability of the complex decreases the pH of the solution needs to be higher.
- It's advisable to use **buffer** solution (a mixture of weak acid and it's conjugate base **OR** vice versa)
- Using buffer will control pH to be nearly constant value

For reactions with M^{2+}

- H_4Y needs NaOH to be added in aqueous solution to make it dissolved.
- Na_2H_2Y is dissolved in aqueous solution
- $M^{2+} + H_4Y \rightleftharpoons MH_2Y + 2H^+$

** If the formation constant is very high and thus equilibrium lies to right

What is the media of reaction?

- ** The reaction should be in a buffer medium (WHY?)

**using basic solution (buffer) will remove H^+ as it's formed so equil. goes to right favoring EDTA-metal complex.

**higher the acidity of the solution leads to reduce the stability of formed complex

Complexometric Indicators:

- Organic dyes form complex with metal weaker than the complex formed by the metal with EDTA
- These indicators possess acid-base indicator properties I
- At the end point indicator will be displaced completely by EDTA, thus the free indicator serves as end point indicator
- Example: Murexide and Eriochrome black T

Indicator that will be used in our exp.:

****Eriochrome black T :**

the color changes from wine red to pure blue

****Murexide :**

the color changes from pink to purple

Objectives:


- Calculate of water hardness (amount of dissolved Ca^{2+} and Mg^{2+} in water
- Calculate conc. of Ca^{2+} and Mg^{2+} in an unknown sample in ppm (mg/L)

a) Standardization of EDTA against $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$ as primary standard


$$M_{\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}} = \frac{\frac{m}{M_w}}{v}, \quad v = 0.1 \text{ L}$$

$$(M \times V)_{\text{EDTA}} = (M \times V)_{\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}}$$

b) Determination of water hardness due to presence of CaCO_3 and MgCO_3



The color changes from wine red to pure blue



The color changes from pink to purple

V1=volume of EDTA corresponding to Ca^{2+} and Mg^{2+}

V2= volume of EDTA corresponding to Ca^{2+} only

V1-V2 =volume of EDTA corresponding to Mg^{2+} because $\text{Mg}^{2+} + \text{OH}^- \longrightarrow \text{Mg}(\text{OH})_2$

$$(M \times V_1) \text{ EDTA} = (M \times V) \text{ Ca}^{2+} \text{ and } \text{Mg}^{2+}$$

$$(M \times V_2) \text{ EDTA} = (M \times V) \text{ Ca}^{2+}$$

$$M \times (V_1 - V_2) \text{ EDTA} = (M \times V) \text{ Mg}^{2+}$$

$$\text{Conc g/L} = M \times \text{Mw}$$

Water hardness In terms of CaCO_3

$$= M \text{ Ca}^{2+} \text{ mol/L} \times 100 \times 1000$$

c) Determination of the concentration of Ca^{2+} and Mg^{2+} in an unknown solution

$$M \times V_2 \text{ EDTA} = (M \times V) \text{Ca}^{2+}$$

$$M \times (V_1 - V_2) \text{ EDTA} = (M \times V) \text{Mg}^{2+}$$

$$\text{Conc g/L} = M \times M_w$$