1 Techniques and Measurements

Name:	Section
Lab. Instructor	Date
Pre-Laboratory Qu	estions
1. A metal sphere weighing 9.48 g is graduated cylinder. The water level density of the metal.	rises to 24.78 mL. Calculate the
2. An empty beaker weighs 32.4257 g liquid is transferred to the beaker. liquid sample was 39.4507 g. Calc liquid.	The total mass of the beaker and culate the density of the unknown
density of the Band :-	70250 - 0. 702 g/ml
3. A term that is easily confused with meant by specific gravity? What are	the units of specific gravity?
volume of the substance equal volume of the Specific growity has	reference substance Lwater
Specific gravity =	al HaO

1 Techniques and Measurements

Name:	Section
Lab. Instructor	

Results and Calculation

A. Determination of the Density of Pure Liquid

Pure water	Trial (I)	Trial (II)
Mass of beaker	32.96 t.0001 g	32.95 ±0.01 g
Mass of beaker + water	42.82±0.01 g	42.83 to.01 g
Mass of water	9.86 ± 0.01 g	9.88 ± 0.01 g
Volume of water	10.0 to.1 mL	100 ± 0.1 mL
Temperature of water	180 ± 0.2 °C	18.0±0.2 °C
Density	0.986 g/mL	0.988 g/mL
Average density	10.9	84 g/mL
Handbook density	0,997	862 g/mL
Unknown liquid		// whole
Unknown Nomber	(A released	eil ni level helen lene
Mass of beaker	32.96 ± 0.01	32.95 ± 0.01 g
Mass of beaker + unknown liquid	45.21 ± 0.01	g 45.22 to.01 g
Mass of unknown liquid	12.25 ± 0.01	g 12.27 to.01 g
Volume of unknown liquid	15.9 ± 0.1 m	L 16.0 ± 0.1 ml
Density of unknown liquid	(0,4404) g/m	nL 0.4668 g/m
	1 /X	

B. Density of Solutions

	Trial (I)	Trial (II)
Mass of beaker	32,96 t 0.01 g	32.95 40.01
Volume of solution	10.0 ± 0.1 mL	10.0 ± 0.1 m
Mass of beaker + solution	43.51 to.01 g	43.51 ± 0.01 g
Mass of solution	10.55 9	10.56
Temperature of solution	19.0 ± 0.2 °C	19.0 ± 0.2 °C
Density of solution	1.055 g/mL	1.056 g/mL

C. Density of Solids

	Trial (I)	Trial (II)
Mass of beaker	33.22 20.019	33.23 + 0:01 g
Mass of beaker + solid pieces	412.173 ± 0.01 g	42.73 £ 0.01 g
Mass of solid pieces	9.51 ± 0.01 9	9.50 ± 0.01 9
Initial water level in the graduated		RELIGIO DE PROPERTORIO
cylinder	25.1 ± 0.11 mL	25.3 ± 0.1 mL
Final water level in the graduated		
cylinder	29,2 ± 0,1 mL	29.1 ± 0.1 mL
Volume of solid pieces	4.1 mL	3.8 mL
Density of the solid substance	2.3/1 g/mL	2.50 g/mL

Questions

1.	What error would be introduced into the density of the metal pellets if you had not shaken the pellets to remove adhering air bubbles? Would the density be too high or too low? The density would be too low because.
	The density of silver is 10.5 g/cm³ and the density of platinum is 21.45 g/cm³. If equal masses of silver and platinum were transferred to equal volumes of water in graduated cylinders, which graduated cylinder would have the highest volume reading? Explain.

2 Formula of a Hydrate

Na	me:
Lat	o. Instructor Section
	Pre-Laboret Date
	Pre-Laboratory Questions
1.	What are alums? Give examples other than potassium alum.
	. Alums are a class of hydrates that consist a
	double salt with the general formula (M'M" (504)2 X H20)
	An example is the potassium alum KAI (SOU)2.12H2O.
2.	What are the hydrates? Give few examples of metal salt hydrates.
	Hydrates are sold compounds that have been
	crystall?zect from aqueous solutions with weakly.
	bounded conter molecules contained in the crystal.
	An example es copper sultate pentahydrate CusouBH2C
3.	Potassium chromic alum has the formula: KCr(SO ₄) ₂ .xH ₂ O. A sample
	of 1.12 g of this alum was heated in a crucible to get a constant mass.
	The mass of the anhydrous salt produced $(KCr(SO_4)_2)$ was 0.64 g.
	o I let the value "y" in the formula of the alum.
	mass of HaO = 1,12q - 0.64g = 0,40g
	moler mass of H20 = 18.002 g/mo!
	n. of H20 moles = 0.48 = 0.0266 mol
	n. of anhydrops sult moles = 0.64 = 2.26 × 103 mol


2 Formula of a Hydrate

Name:	Section
Lab. Instructor	Date

Results and Calculations

A. Potassium Alum:

Mass of empty crucible (m ₁)	16.48 ±0.01	g
Mass of crucible and the alum (m ₂)	17.78 ±0.01	g
Mass of crucible and anhydrous salt (m ₃)	14.35 ± 0.01	g
Mass of alum (m ₂ -m ₁)	1.00	g
Mass of anhydrous salt (m ₃ -m ₁)	0.57	g
Mass of water lost upon heating (m2-m3)	0.43	g
Number of moles of water lost upon heating	18.016 = 0.021	11101
Number of moles of anhydrous salt (KAI(SO ₄) ₂)	0.5H 2.21 * 258.22	
Percentage of water of crystallization, by mass alam	1.00 \$ 100 =	43%
The value "x" in the formula, (number of moles of	0.024	10.8
water of crystallization / number of moles of	2.21*103	11
anhydrous salt)		



QUESTIONS

1 What is the offers
1. What is the effect on the calculated value of "x" if the dehydration of the alum is not completeThevalueoF
2. A student heated 1.16 g of hydrated sodium sulfate in a crucible to get
0.51 g of anhydrous salt. What is the formula of the anhydrous salt? (Show your work)
mass. of water = 0.65 g
n. moles of water = 0.036 ma)
3.59 × 10 ³
Formula: Na2804.10 H20
10

3 The Empirical Formula of an Oxide

Name: Section
Lab. Instructor Date
Pre-Laboratory Questions
1. What is the mole and molar mass? What are the molar masses of magnesium and atomic oxygen? Mole 9s a work of measurement for the amount of substance or cheroscal amount. Molar mass 9s a physical property characteristic or
2. When 0.192 g of phosphorus is burned, 0.341 g of a white oxide is obtained. (a) Write a balanced chemical equation for the reaction of phosphorus with molecular oxygen based on this empirical formula. (b) Determine the empirical formula of the oxide.
a) 4? 302 > 29203 b) n. phospharus = 0.192 = 6.19 * 103 mol 30.94
$0.02 = 0.149.9.$ $0.02 = 9.31 * 10^3 \text{ mol.}$ $0.02 = 9.31 * 10^3 \text{ mol.}$ $0.03 = 0.149.9.$ $0.03 = 0.149.9.$ $0.03 = 0.149.9.$ $0.03 = 0.149.9.$ $0.03 = 0.149.9.$ $0.03 = 0.149.9.$ $0.03 = 0.149.9.$ $0.03 = 0.149.9.$ $0.03 = 0.149.9.$ $0.03 = 0.149.9.$ $0.03 = 0.149.9.$ $0.03 = 0.149.9.$ $0.03 = 0.149.9.$ $0.03 = 0.149.9.$ $0.03 = 0.149.9.$ $0.03 = 0.149.9.$ $0.03 = 0.149.$
$(P,O_{15})*2 \Rightarrow P_2O_3$ emprical formula

3 The Empirical Formula of an Oxide

Name:	Section
lab instructor	Date

Results and Calculations

Mass of empty crucible (after first heating)	15.64 g	
Mass of empty crucible (after second heating)	15.63 g	
Final mass of empty crucible	15.63 g	
Mass of crucible and Mg	15.83 g	
Mass of Mg	0.20 g	
Moles of Mg (n1)	0.20 8.23 ± 10 mol	3
Final mass of crucible and Mg-oxide	15,94 9	
Mass of Mg-oxide produced	0.31	
Mass of oxygen gained	el 6.71 g	
Moles of oxygen atoms (n2)	6.87 * 103 mol	
Formula of magnesium oxide (Mg _{n1} O _{n2})	My 8,73 × 103 6	87.1
Empirical formula of magnesium oxide	MgO	
Mass percent of Mg in the oxide (x ₁) (experimentally)	64.5 %	
Mass percent of Mg in the oxide (x2) (calculated for MgO)	60.3 6%	
Percentage error = $[x_2-x_1 /x_2] \times 100\%$	6.9 =7 %	

QUESTIONS

1.	If water had not been added to your initial product, what error in the determined percentage of magnesium would have resulted (that is, if part of the product has been magnesium nitride)? Explain.
	The percent would Properse because the mass of Mg would Properse and mass
	of MgO would decrease
	es al
2.	If large amount of magnesium oxide had been lost during the heating
	of the crucible, would this have made the calculated %Mg in the product too high or too low? Explain.
	Too high because the mass of MgO
om.	would decrease

4 Limiting Reactant

Lab. Instructor	Section
	Date
Pre-Laboratory Q	uestions
1. Na ₃ PO ₄ .12H ₂ O and BaCl ₂ .2H ₂ O are exampled hydrates. Calculate the molar masses of Nand anhydrous forms.	amples of compounds known as la_3PO_4 and $BaCl_2$ in the hydrated la_3PO_4 and $BaCl_2$ in the hydrated la_3PO_4 and l
n. Bacla = 13 x 103 mol	n. NazPoy . 4.4 * 10 mg1
m. Baz (POw) 2 produced from	Bacl2 1-
1.3 * 103 mol Bacl2 * 1 m	ol (Bas(POub) = 43 * 10 mol
+ n Buz (PO4) 2 produced From	Nazpoy:
4.4 * 163 mol Daz Pou * 1 mo	of Bact2 = 2.2 * 103
80 Baclz Ps the 55 19mg	Himy reagent
, m. Bag (Pay)2 = 601.84 * 4.3	3 *10 = 0.22 9

4 Limiting Reactant

Results and Calculations

A. Precipitation of Ba₃(PO₄)₂:

Unknown Number:

Mass of salt mixture (m ₁)		
Mass of filter paper (m ₂)	0.40	9
	0.68	g
Mass of filter paper and Ba ₃ (PO ₄) ₂ (m ₃)	0.84	g

B. Determination of the Limiting Reactant:

Limiting reactant in salt mixture is3.

Excess reactant in salt mixture is

Ba.Cl, 2H20 Na3Po.4.12H20

Mass of Ba ₃ (PO ₄) ₂ precipitated (m ₃ -m ₂)	0.29	g
Number of moles of Ba ₃ (PO ₄) ₂ precipitated (n₁)	4.84 10	mol

(a) If the limiting reactant is BaCl₂. 2H₂O:

Number of moles of BaCl ₂ .2H ₂ O reacted (n ₂)	1.41 × 10 ³ mol
Number of moles of Na ₃ PO ₄ .2H ₂ O reacted (n ₃)	9.62×104 mol
Mass of BaCl ₂ .2H ₂ O reacted (m ₄) N2 × MW	o. 34 g
Mass of Na₃PO₄ reacted (m₅) √3⊀ Mω	0,160.36 9
Mass of excess Na ₃ PO ₄ [m ₁ -(m ₄ +m ₅)]	0.19 Zeno 9
Mass percentage of BaCl ₂ .2H ₂ O	493 %

my x100%

QUESTIONS

1. Calculate the mass of Ba ₃ (PO ₄) ₂ produced from the reaction of 0.78 g Na ₃ PO ₄ .12H ₂ O with excess BaCl ₂ .2H ₂ O. What is the purpose of heating the mixture in step 3 for 20 minutes? In order for the reaction to (coagulation) mm. Na ₃ PO ₄ 12H ₂ O = 380.2 g/me.
D. DagPoy 12H20 = 0.78 = 2.1 * 103 mal.
380.2 * 10 mal.
$2.1 \times 10^{-3} \text{ mol} \times 1 \text{ mol} = 1.02 \times 10^{3} \text{ mol} \text{ Bas (POL)}_{2}$
m Baz (PO4)2 = 601.84 * 1.02 * 103 = 0.61 a
2. What is the purpose of washing the precipitate with hot water in step 4?
How would the reported percentage of the excess reactant be affected
if the precipitate was not washed in this step?
To collect all the preceptate. The percentage
of the strategy readent could francise
causing the excess reagent to decrease
(the 19miting reactant cost be a mixture of
both LB and EB)

5 Determination of Acetic Acid in Vinegar

Lab. Instructor	Section Date
Pre-Laboratory C	Questions
1. Why was the standard NaOH se	olution not prepared by calculating ed for 100 mL of solution, weighing it xactly 100 mL of total volume?
which means that 9ts consent	trateon changes with the
and 9t may react with	
2. Why not simply rinse the buret wis solution to be used in it? This assures that the solution buret will not be different water adhering to the	than to be used an the
3. Why does the volume of water ac	
Because the basis of this	
one reactant leacts complete	tely with another reactant
with no excess remaining.	
of each reactant is the	
	61

5 Determination of Acetic Acid in Vinegar

Name:	Section
Lab. Instructor	Date

Results and Calculations

A. Standardization of NaOH Solution

Manual Control	Trial (I)		Trial (II)	
Mass of flask	122.08	g	122.08	g
Mass of Flask + KHP	122,28	g	122.29	g
Mass of KHP	0.20	g	0.21	g
Molar mass of KHP	204.22	g/mol	204.22	g/mol
Moles of KHP	98 * 104	mol	(C1.00* 103	mol
Initial buret reading	20.00	mL	18.16	mL
Final buret reading	24.00	mL	124.10	mL
Volume of NaOH	M.00	mL	6.00	mL
Moles of NaOH	9.8 * 104	mol	1.0 * 103	mol
Molarity of NaOH	0.14	M	0.14	М
Average Molarity of NaOH		0.16		M

B. Mass percent of acetic acid in vinegar. (15)

		Fare to	
	Trial ((1)	Trial (II)
Volume of vinegar	10.00	mL	10.00
Initial buret reading	22.00	mL	20.80
Final buret reading	36.60	mL	35:00
Volume of NaOH used	141.60	mL	14.20
Average molarity of NaOH, From P	art A:		0.16
Moles of NaOH used n=cv	2.3	mol	2.3 m
Moles of CH ₃ COOH in vinegar Reacted with the NaOH	2.3	mol	2.3 m
Molarity of CH ₃ COOH in vinegar	0.923	M	0.023
Average molarity of acetic acid in vi		eve /	0.23
Molar Mass of acetic acid, CH ₃ COO	Hm: M&n	/60	0.05% g/m
Mass of CH ₃ COOH per liter of vineg	ar	14 131	3 * 103
Mass percent acetic acid in vinegar			
assume that vinegar has a density	of 1.00 g/mL)	0.1	4
$d = \frac{m}{v}$			
m= 1.00 + 10.00			
= 10.00 a			
1			

Generated by CamScanner from intsig.com

QU	-0		0		(04
		6 B I		V. 1	
E4 1 E 1					

Ireal - [] + 100 %

Calculate the percent error would have been in a titration that used
32.75 mL of a solution if a bubble with a volume of 0.25 mL had
heen swept out of the tip during the titration.

0.25 × 100% = 0.46% 32.75 32.50

The label on the vinegar bottle used in this experiment claims that the vinegar contains 3 % acetic acid by weight. Use your results and a density of 1.0 g/mL to investigate this claim.

3.0/0 = Conc = 60.05? = 1000/0

concentration = a 499 / 0.5 M

6 The Neutralizing Capacity of Antacid Tablets

Name: Lab. Instructor	Section Date
Pre-Laboratory Ques	
1. Write a balanced equation for the reaction.	on of HCI and NaHCO ₃ .
2. How many moles of HCI are needed to	
	recieve
tablet? Because Strong bases as action bases	go to completion
and so this case we don	it count the stomach
to become neutralized or. Postead we need 9t to normal acrd9ty.	get back to Pt's

6 The Neutralizing Capacity of Antacid Tablets

Name:	Section
. L Instructor	Date

Results and Calculations

\$10/436 In \$0.04 lanis) o.e	Trial (I)	ntered	Trial (II)	
Mass of antacid sample	0.20	9	0.21	9
Volume of 0.150 M HCI solution	25.00	mL	25.00	mL
Moles of HCl (used to dissolve antacid)	3.75 * 10	mol	3.75 * 103	mol
Initial buret reading	0.00	mL	0.00	mL
Final buret reading	16.20	mL	16.40	mL
Volume of NaOH added	16.20	mL	16.40	mL
Moles of NaOH (used to titrate the excess acid)	2.43 * 10	³ mol	2.46 × 103	mol
Moles of excess HCI	2.43 * 10	mol	2.46 * 163	mol
Moles of HCI(needed to neutralize the antacid tablet)	1.32 * 103	mol	1.20 * 103	mol
Neutralizing capacity of antacid	6.60 × 103 n		6.141 × 103 mc	HCI/g
Average (mol HCI/g antacid)	6.8H	* 103	mol HCl/g	antacid

QUESTIONS

1. Assume the concentration of stomach acid is 0.10M and the neutralizing capacity of the acid is 0.013 mol/g, calculate the volume of stomach acid needed to neutralize 1.0 g antacid tablet.
2. A 0.333 g-sample of antacid was dissolved in 40.00 mL of 0.135 M
HCI solution, then back-titrated to the end-point with 9.28 mL of a
0.0203 M NaOH solution.
a) Calculate number of moles of acid in the original 40.00 mL of HCI?
n. moles: 40 * 103 * 0.135
= 5.46 × 103 mol
b) How many moles of base were used in the back titration of excess/
ask & 1 & 1 & 1 & 1 & 1 & 1 & 1 & 1 & 1 &
- 1 & 8.9 . * . ! S
c) How many moles of excess HCI !?
1.88 ×10 mol
W.9.
d) How
d) How many moles of HCI reacted with the antacid sample?
$= 5.21. \pm 10^3 \text{ mol}$

Date
Pre-Laboratory Questions
1. Explain why cooked fruits and vegetables have lower vitamin C
Cooking leaches water solluble returnin C. From the vegetab
and high temperatures accelerate the degradation by air
What are the oxidizing agents in this experiment?
KI, Is 703
What will be the color change of the starch indicator that indicates the end of the titration in this experiment?
Vitamin C is an acid (ascorbic acid) and a reducing agent. Which
property is utilized for this analysis in this experiment?
redocing eigent
If 39.50 ml of a vegetable juice contains 45% of the recommended
daily allowance of vitamin C (equal to 60 mg). How many milliliters of
the vegetable juice will provide 100% of the recommended daily
allowance?
39.50 ml 0.45
V = 84.40 mL

7 Vitamin C Analysis

Name:	Section
1 - structor	Date

Results and Calculations

	-	
	Trial (I)	Trial (II)
Mass of sample	0.10 g	9
Volume of 0.010 M KIO ₃ added	50.0 mL	mL
Moles of IO ₃ added	5 * 10 4 mo	mol
Moles of la generated, total n & 3	1.5 *103 mo	mol
Initial buret reading	11.00 ml	mL mL
Final buret reading	24.00 ml	mL mL
Volume of 0.10 M Na ₂ S ₂ O ₃ added	16.0 ml	ML mL
Molar concentration of Na ₂ S ₂ O ₃	0-1 mol/l	_ mol/L
Moles of S ₂ O ₃ ²⁻ added	1.6 × 10 mc	mol mol
Moles of l_3 reduced by $S_2O_3^{2-}$ $(\frac{\Omega}{2})$	8 * 10 mc	mol mol
Moles of I ₃ reduced by C ₆ H ₈ O ₆	74104 ma	ol mol
Moles of C ₆ H ₈ O ₆ in the sample	14 10 mg	ol mol
Mass of C ₆ H ₈ O ₆ in the sample	0.12	g g
Mass percent of C ₆ H ₈ O ₆ in the sample	123	%
Average percent of C ₆ H ₈ O ₆ in the sample	13	

AI			0	ALC: U	-		
	B 2	_		E 86			LN
2						LBC.	

1. A 25.0 mL volume of 0.010 M KIO ₃ , containing an excess of KI, is added to a 0.346 g sample of a lemon solution containing vitamin c. The red brown solution, caused by the presence of excess I ₃ , is titrated to a colorless starch end point with 7.70 mL of 0.100 M Na ₂ S ₂ O ₃ .
a. Calculate the moles of l ₃ that are generated from KIO ₃
titration?
$\begin{array}{cccccccccccccccccccccccccccccccccccc$
lemon sample?
ask
d. Calculate the moles and grams of vitamin C in the sample.
mass.=
e. Calculate the percent (by mass) of vitamin C in the lemon sample. 2. If the blue color does not appear when the starch solution is added during the titration, should you continue titrating or discard the sample? Explain.
Discard because this means that the end-point was passed 88

Generated by CamScanner from intsig.com

Name:	Section
Lab. Instructor	Date
Pre-Laboratory Quest	tions
1. Write balanced equation(s) for the reaction of sodium thiosulfate. 6.5.2.03. L. Cr. 20.7. L.	
2. In today's chemical analysis of bleach oxidized by the hypochlorite ion? Lockete Formulas of sodium thiosulf	
4. Does sodium thiosulfate serve as an reducing agent? What does sodium thios	sulfate reduce?

Name:	Section
teuctor	Date

Results and Calculations

A. Standardization of thiosulfate solution

	Trial (I)		Trial (II)	
Volume of 0.010 M K ₂ Cr ₂ O ₇ solution	25.0	mL	25.0	mL
Initial buret reading	0.00	mL	15.0	mL
Final buret reading	15.0	mL	30.2	mL
Volume of Na ₂ S ₂ O ₃ solution	15.0	mL	15.2	mL
Number of moles of K ₂ Cr ₂ O ₇	2.8 4104	mol	cieve 10	mol
Number of moles of Na ₂ S ₂ O ₃ $\cap *$ 6	1.5 * 103	mol	1.5 4 103	mol
Molarity of Na ₂ S ₂ O ₃	0.1	M	6.1	M
Average molarity of Na ₂ S ₂ O ₃		0.11	6	N

B. Analysis of bleach solution

	Trial (I)		Trial (II)	
Dilution Factor	20		20	
Volume of bleach solution	25.0	mL	25.0	ml
Initial buret reading	21.0	mL	31.2	mL
Final buret reading	31.2	mL	41.0	mL
Volume of Na ₂ S ₂ O ₃ solution	10.2	mL	MOTON	mL
Number of moles of Na ₂ S ₂ O ₃ oxidized	1.02 * 103	mol	1.01 * 153	mol
Number of moles of CIO reduced (2)	8.10 = 10 4	mol	5.05 * 10	⁴ mol
Molarity of diluted bleach solution	0.0209	M	0.0202	M
Molarity of original bleach solution	9.408	M	0.404	M
Average molarity of original leach solution	& re	cie	ve /	M
Mass% of NaClO(assume ensity of bleach solution = .00g/mL)		3.02		%
- M orrgral = Malil *				

QUESTIONS

1. A 10.0 mL bleach sample is diluted to 100 mL in a volumetric flask. A 25.0 mL of this solution is analyzed according to the procedure in this
experiment. If 12.5 mL of 0.30 M Na ₂ S ₂ O ₃ is needed to reach the end
point, calculate the mass percent of NaClO in the original sample?
(Assume the density of bleach solution is 1.084 g/mL).
n. Neu25203 + 3.75 * 103 mel
0. $C_{10} = (3.75 * 10^3) = 1.88 * 10^3$
4 = 0.018 mol/L
ach so of recieve
as-/ mass - 25/22 % recieve
2. A bleach solution that is 6.58% NaClO (density= 1.10 g/mL) was diluted
to 0.056 M concentration Calculate dilution factor.
% (F8 - M + 74.44 100 %
% 6.58

Hi. F. Q. 97
D.F = H1 = 17.32
MF

Molar Mass of a Volatile Liquid

Name:	Section
Lab. Instructor	Date
Pre-Laboratory Quest	ions
1. Dumas method assumes the gas behave when gases behave ideally? At high temperatures and	
a) What is the molar mass of the volatile	The mass of vapor is the start of the start
b) What is the density of vapor at STP? PV = n.RT PV = m.RT MP = m.RT	h= 100m

9

Molar Mass of a Volatile Liquid

Name:	Section
Lab. Instructor	Date

Results and Calculations

	Trial (I)		Trial (II)
Boiling point of water	94.6	°C	°C
Atmospheric pressure	682	mmHg	mmHg
Mass of empty flask	M8.53	g	g
Mass of flask and condensed vapor	1918	g	9
Mass of condensed vapor	50.65	g	9
Volume of flask	206	mL	mL
Boiling point of water	E340.8	reg	eve / ĸ
Atmospheric pressure	0.89M	atm	atm
Volume of flask	0.195	L	L
Gas constant (R)			0.0821 L.atm /mol.K
Molar mass of unknown	113	g/mol	g/mol
Average molar mass			g/mol

QUESTIONS

1	If the outside of the flask is not dried after vaporizing the liquid, will the calculated molar mass be too high or too low? Explain.
	To blah a because the mass of the
	Flusk would Increase due to the water on the walls of the Flask
	Consider the following experimental data: Mass of condensed vapor =
2.	0.395 g at 96 °C and 755 mm Hg occupies 137 mL. What is the molar mass of the liquid?
	= 0.398(0.0821)(953242)
	= 0.395 (0.0821) (95.4243) (0.134) (455 1450)
	(0.134) (455 1450)

10	ary aria ricco o Law
Name:	. Section
Lab. Instructor	Date
Pre-Laboratory C	luestions
2. In a calorimeter calibration experate at 55.2 °C is added to a calorin 23.5 °C. After stirring and waiting final temperature reached is constant.	riment, a sampe of 51.682 g of water neter containing 50.220 g of water at ing for the system to equilibrate, the 36.6 °C. Calculate the calorimeter
3. Give chemical equations for the	reactions that will be a
	Nacl. 1. H20. 1. UH3

10 Thermochemistry and Hess's Law

Name:	Section
Lab. Instructor	Date

Results and Calculations

A. Determination of a Calorimeter Constant.

	Trial (I)		Tri	ial (II)
Mass (or volume) of cold water	75	g		g
Temperature of cold water	25.0	°C		°C
Mass (or volume) of hot water	45	9		g
Temperature of hot water	70	°C		°C
Final temperature reached	46.0	°C		°C
Increase of cold water temperature (ΔT)	21.0	°C		°C
Decrease of hot water temperature (ΔT)	regie	vêc	11	°C
Heat lost by hot water	4531.2			J
Heat gained by cold water	6590	1		J
Heat gained by calorimeter water	941.2/	/	J	J
Calorimeter constant (C)	44.8	J/°(1/°C
Average value of calorimeter constant	1			1100

- qwarm = qcold + qcal

B. Heat of Acid/Base Reactions

1. HCI/NaOH

	Trial (I)		T.
Volume of 1.0 M NaOH used	50	mL	Trial (II)
Initial temperature of NaOH	22.0	°C	
Volume of 1.0 M HCl used	50	mL	0
Initial temperature of HCI	22.0	°C	
Final temperature reached	24.5	°C	
Total volume of mixture	100	mL	0
Total mass of mixture	100	g	n
Temperature change, Δt	5.5	°C	0
Average value of the Calorimeter constant (C)	44.8		J/0
Heat gained by the solution*	2239	J	
leat gained by the calorimeter	0 246	J	
leat of the reaction	- 2485	J	/
loles of NaOH reacted		not	
oles of HCI reacted	1	nol	mo
oles of water produced	7	nol	
2485	/		kJ/mol H₂ ⁰
erage value of ΔH	7		kJ/mol H ₂ O
erature value of ΔH	49.4		- 55.9 kJ/mol H ₂ 0

^{*} Assume density of solution = 1.00 m/m

2. HCI/NH3			
Whime of 1 O M NULL	Trial (I)		Trial (II)
Volume of 1.0 M NH ₃ used	50.0	ml	Tital (II)
Initial temperature of NH ₃			mL
Volume of 1.0 M HCl used	22.0	°C	°C
Initial temperature of HCI	50.0	mL	mL
	22.0	°C	°C
Final temperature reached	24.0	°C	
Total volume of mixture	2.1.0		°C
Total mass of mixture	100	mL	mL
	100	g	g
Temperature change, ∆t	5.0	°C	°C
Average value of the Calorimeter constant (C)	44.8		J/°C
Heat gained by the solution*	2035	J	J
Heat gained by the calorimeter	2241	J	J
Heat of the reaction	- 2259	J	1350 F J
Moles of NH ₃ reacted	0.05	mol	mol
Moles of HCI reacted	0.05	mol	mo
Moles of NH ₄ CI produced	0.05	mol	mo
Н	45.2 kJ	/mol	kJ/mo
verage value of ΔH	1	45.2	kJ/mo
H for the reaction:	4.5	ų.	k
aOH + NH ₄ CI → NaCI + NH ₃ + H ₂ O	1.0		

^{*} Assume density of solution = 1.00 g/mL and specific heat of solution = 4.07 J/g. °C.

Questions

1.	What effect on the calculated calorimeter constant would be	Ob
	if the calorimeter cup were made of conducting material	Observed
	metal) rather than plastic foam?	(such as

It would movedse

2. Why is water typically used as the heat-absorbing liquid in calorimeter?

Because water has a very high specific heat capacity than other liquids

3. Given that:

$$H_2O_{2(\ell)} \rightarrow H_2O_{(\ell)} + \frac{1}{2}O_{2(g)}$$
 $\Delta H = -96.0 \text{ kJ/mol}$ $2H_{2(g)} + O_{2(g)} \rightarrow 2H_2O_{(\ell)}$ $\Delta H = -571.6 \text{ kJ/mol}$) %2

Use Hess's law to calculate ΔH for the reaction:

$$H_{2(g)} + O_{2(g)} \rightarrow H_2O_{2(g)}$$

H2 1/202 H202 AH = 96.0 KJ [mol H2 + 02 L7 H202 AH = 1888 KJ [mol

11 Molar Mass from Freezing Point Depression

Name: Section
Lab. Instructor Date
Pre-Laboratory Questions
1. Define what we mean by the colligative property.
collegative property in are properties of solution.
that depend on the
that depend on the number of molecules in a
gover volume of solvent and not on the properties
- Polentity (59ze/mass) of the molecules
2. A 0.36-g sample of an unknown substance was dissolved in 30.0
mL of cyclohexane. The density of cyclohexane is 0.779 g/mL.
The freezing-point depression was 2.50 °C. Calculate the molar
mass of the unknown substance. [k _f (cyclohexane) = 20.0 °C/m]
O ATE = KEM mass
2.5 = 20 m n.s. 0.105. * 23.3.K. * 16
m = 0.125 mayky = 2.921 *153 mol

0 d= m
m = 0.449 * 30 M = 0.36
= 23.349
=23.34 × 163 Kg = 123.2
25.51 2 10 g/nels

11 Molar Mass from Freezing Point Depression

Name:	Section
I ah Instructor	Date

Results and Calculations

A. Freezing point of cyclohexane

	Trial (I)	Trial (II)	
Mass of empty small tube	29.84 g	29.84	g
Freezing point of cyclohexane	7.0 t 0.2 °C	7.0 to.2	°C
Average value of freezing point	1	1.0	°C

Hw = KI * msoilute

B. Freezing point of solution

	Data	
Mass of empty small tube	29.84	g
Mass of solute	0.21	g
Mass of empty small tube + cyclohexane	34.49	g
Mass of cyclohexane	7.65	g
First freezing point	4.0	°C
(ATf)1 (M.O-4.0)	3.0	°C
Molar mass of solute	183	g/mol
Mass of second sample of solute	0.15	g
Total mass of solute	0.36	°C
Second freezing point	3.0	10
$(\Delta T_f)_2$	4.0	
Molar mass of solute	134	
Avon		

Generated by CamScanner from intsig.com

QUESTIONS 1. What would be the effect of each of the following on the calculated molar mass of the solute?
a) Some cyclohexane evaporated while the freezing point of pure cyclohexane was being measured. no effect
b) Some cyclohexane evaporated after the solute was added would decrease
c) A foreign solute was already present in the cyclohexane.
that is 1.5 °C too low at all temperature.